



MP Biomedicals, LLC

29525 Fountain Parkway
Solon, Ohio 44139

Telephone: 440/337-1200
Toll Free: 800/854-0530
Fax: 440/337-1180
mailto: biotech@mpbio.com
web: <http://www.mpbio.com>

TECHNICAL INFORMATION

ENZYME SYSTEMS PRODUCTS

a division of



Technical Data Sheet

Method for Assay of Caspase 2 with Ac-Val-Asp-Val-Ala-Asp-AFC

Materials:

100 mM HEPES, pH 7.5, 20% (v/v) glycerol, 5 mM DTT, 0.5 mM EDTA

– Buffer

20 mM stock solution of Ac-Val-Asp-Val-Ala-Asp-AFC in DMSO

– Substrate

Cell lysate or purified enzyme solution (~15 nanograms enzyme)

– Enzyme

80 μ M free AFC (Catalog # T07) in DMSO

– Fluorescence Standard

Method:

- Add 10 μ l of enzyme to 490 μ l buffer. Mix. Incubate at 30° C for 30 minutes.
- With fluorometer adjusted to 400nm excitation and 505nm emission, add 20 μ l of substrate to enzyme solution.
- Record increase in fluorescence from T_0 to T_{end} where fluorescence generated at T_{end} are significantly different from those at T_0 .
- Record fluorescence units generated by 10, 20, and 30 μ l free AFC 490, 480, and 470 μ l buffer solution, respectively.
- Graph fluorescence units vs. micromole AFC. Use slope to convert fluorescence units generated by enzyme to activity.

Storage:

Desiccate AFC142 in solid form at room temperature. Store DMSO/DMF solutions at -20° C. Material is stable for at least one year, if stored as recommended.