

# HIV BLOT 2.2

## WESTERN BLOT ASSAY

### Instructions For Use

**0123**  
 REVISION DATE 2016-05  
 MAE0011-ENG-5

(18 tests kit) : 11030-018  
 (36 tests kit) : 11030-036

#### NAME AND INTENDED USE

The **MP Diagnostics HIV BLOT 2.2** is a qualitative enzyme immunoassay for the *in vitro* detection of antibodies to human immunodeficiency virus type 1 (HIV-1) and type 2 (HIV-2) in human serum or plasma. It is intended for use as a more specific supplemental test on human serum or plasma specimens found repeatedly reactive using screening procedures such as the Enzyme-Linked Immunosorbent Assay (ELISA).

#### INTRODUCTION

Screening tests are widely available for detecting antibodies to both HIV-1 and HIV-2, the etiologic agents of the Acquired Immunodeficiency Syndrome (AIDS). Such tests can be extremely sensitive but have a potential for being less specific, leading to false positive interpretations. Independent supplemental tests of high specificity are therefore necessary to further confirm the presence of antibodies to HIV-1 and/or HIV-2.

The **MP Diagnostics HIV BLOT 2.2** kit is intended for use as a more specific supplemental test on human serum or plasma specimens found repeatedly reactive using ELISA. The separated specific HIV-1 viral antigens incorporated onto the strips via electrophoretic and electrotransblot procedures, combined with a specific HIV-2 synthetic peptide on the same strip allow for further delineation of the antibody responses to specific viral proteins. Each strip also includes an internal sample addition control to minimize the risk of false negatives due to operational errors and to ensure the addition of samples.

#### REAGENTS

- Loosen cap of sample container.
- Heat-inactivate sample at 56°C for 30 minutes in a water bath.
- Allow sample to cool down before retightening cap.
- Sample can be stored frozen until analysis.

#### REPEATED FREEZE-THAWING OF SAMPLE IS NOT RECOMMENDED.

#### ADDITIONAL MATERIALS REQUIRED BUT NOT PROVIDED

- Deionized or distilled water
- Disposable gloves
- Rocking platform (designed with a rocking speed range of 12 to 16 cycles per minute, and which moves through a 5° to 10° tilt to wash membranes evenly)
- Pipettors and tips of appropriate volume
- Aspirator with sodium hypochlorite trap
- 56°C water bath (optional)
- Sodium hypochlorite for decontamination

#### PREPARATION OF REAGENTS

- DILUTED WASH BUFFER**
  - DILUTED WASH BUFFER should be **prepared fresh prior to use**.
  - Dilute 1 volume of WASH BUFFER CONCENTRATE (20x) with 19 volumes of reagent grade water. Mix well.
- BLOTTING BUFFER**
  - BLOTTING BUFFER should be **prepared fresh prior to use**.
  - Dilute 1 volume of STOCK BUFFER CONCENTRATE (10x) with 9 volumes of reagent grade water. Mix well.
  - Add 1 g of BLOTTING POWDER to every 20 ml of the diluted STOCK BUFFER prepared in step 2(b) above. Stir to ensure powder dissolves completely.
  - Stir again before dispensing.
- WORKING CONJUGATE SOLUTION**

Note : Prepare solution in polypropylene container / beaker.

  - WORKING CONJUGATE SOLUTION should be **prepared fresh prior to use**.
  - For **RAPID ASSAY PROTOCOL**, prepare WORKING CONJUGATE SOLUTION by diluting CONJUGATE at 1:500 into BLOTTING BUFFER, for example, 10 µl CONJUGATE to 5ml BLOTTING BUFFER.
  - For **OVERNIGHT ASSAY PROTOCOL**, prepare WORKING CONJUGATE SOLUTION by diluting CONJUGATE at 1:1000 into BLOTTING BUFFER, for example, 5 µl CONJUGATE to 5ml BLOTTING BUFFER.
- SUBSTRATE SOLUTION (ready to use)**
  - Dispense directly the required volume from the bottle. Use a clean pipette. Cap tightly after use.

#### DESCRIPTION OF SYMBOLS USED

	Use by <i>Synonym for this :</i> Expiry Date		In vitro diagnostic medical device
	Batch Code <i>Synonyms for this are :</i> Lot Number Batch Number		Catalogue Number <i>Synonyms for this:</i> Reference number Re-order number
	Temperature Limitation		Cautions
	Manufacturer		Authorised Representative in the European Community
	Contains sufficient for <n> strips		Consult Instructions for Use
	Contents		

#### CHEMICAL & BIOLOGICAL PRINCIPLES OF THE PROCEDURE

The nitrocellulose strips are incorporated with separated, bound antigenic proteins from partially purified inactivated HIV-1 using electrophoretic blotting, plus a specific HIV-2 synthetic peptide on the same strips. Individual nitrocellulose strips are incubated with diluted serum or plasma and controls. Specific antibodies to HIV-1 and HIV-2 if present in the specimens will bind to the HIV-1 proteins and HIV-2 peptide on the strips. The strips are washed to remove unbound materials. Antibodies that bind specifically to HIV proteins can be visualized using a series of reactions with goat anti-human IgG conjugated with alkaline phosphatase and the substrate BCIP/NBT. This method has the sensitivity to detect marginal amounts of HIV specific antibodies in serum or plasma.

#### KIT COMPONENTS

Component Description Provided	Quantity
<b>NITROCELLULOSE STRIPS</b> Incorporated with HIV-1 viral lysate, a specific HIV-2 envelope peptide and a serum addition control band. Keep dry and away from light.	Available in 18 or 36 strips
<b>NON-REACTIVE CONTROL</b> Inactivated normal human serum non-reactive for Hepatitis B surface antigen (HBsAg), antibodies to HIV-1/2, and anti-HCV. Contains sodium azide and thimerosal as preservatives.	1 vial (80 µl)

#### STRONG REACTIVE CONTROL

Inactivated human serum with high titered antibodies to HIV-1 and HIV-2 and non-reactive for HBsAg & anti-HCV. Contains sodium azide and thimerosal as preservatives.

#### WEAK REACTIVE CONTROL

Inactivated human serum with low titered antibodies to HIV-1 ONLY and non-reactive for HBsAg, anti-HIV-2 and anti-HCV. Contains sodium azide and thimerosal as preservatives.

#### STOCK BUFFER CONCENTRATE (10x)

This buffer with heat inactivated normal goat serum. Contains thimerosal as preservative.

#### WASH BUFFER CONCENTRATE (20x)

This with Tween-20. Contains thimerosal as preservative.

#### CONJUGATE

Goat anti-human IgG conjugated with alkaline phosphatase. Contains sodium azide as preservative.

#### SUBSTRATE

Solution of 5-bromo-4-chloro-3-indolyl-phosphate (BCIP) and nitroblue tetrazolium (NBT).

#### BLOTTING POWDER

Non-fat dry milk

#### INCUBATION TRAY\*

Instructions For Use

Forceps

Note : Volume of reagents provided are sufficient for 4 runs.

\* Incubation trays provided but packed separately from the kit.

#### WARNINGS AND PRECAUTIONS

- For *in vitro* diagnostic use only.
- For Professional use only.
- Please refer to the product labelling for information on potentially hazardous components.

#### HEALTH AND SAFETY INFORMATION

**CAUTION:** This kit contains materials of human origin. No test method can offer complete assurance that human blood products will not transmit infection. **HANDLE ASSAY SPECIMENS, STRONG REACTIVE, WEAK REACTIVE AND NON-REACTIVE CONTROLS AS POTENTIALLY INFECTIOUS AGENTS.** It is recommended that the components and test specimens be handled using good laboratory working practices. They should be disposed of in accordance with established safety procedures.

The **Strong Reactive Control, Weak Reactive Control** and **Non-Reactive Control** contain Thimerosal and Sodium azide while **Stock Buffer Concentrate** and **Wash Buffer Concentrate** contain Thimerosal and Conjugate contains Sodium azide. Sodium Azide can react with copper and lead used in some plumbing systems to form explosive salts. The quantities used in this kit are small, nevertheless when disposing of azide-containing materials they should be flushed away with relatively large quantities of water to prevent metal azide buildup in plumbing systems.

Pursuant to EC regulation 1272/2008 (CLP), hazardous components are classified and labelled as follows:

<b>Component:</b>	<b>Nitrocellulose strips</b>
<b>Signal Word:</b>	<b>Danger</b>
<b>Pictogram:</b>	
<b>Hazard Statements:</b>	H228 Flammable solid
<b>Precautionary Statements:</b>	P210 Keep away from heat/sparks/open flames/hot surfaces. – No smoking. P280 Wear protective gloves/protective clothing/eye protection/face protection.
<b>Supplemental Statements:</b>	EUH210 Safety Data Sheet is available on request
<b>Contains:</b>	100% Nitrocellulose

<b>Component:</b>	<b>STOCK BUFFER CONCENTRATE (10x)</b> <b>WASH BUFFER CONCENTRATE (20x)</b>
<b>Signal Word:</b>	Warning
<b>Pictogram:</b>	
<b>Hazard Statements:</b>	H373 May cause damage to organs through prolonged or repeated exposure

samples are preferred. Lipemic, icteric or contaminated (particulate) samples should be filtered (0.45µm) or centrifuged before testing.

Samples can be inactivated but this is not a requirement for optimal test performance.

Inactivate as follows:

- Loosen cap of sample container.
- Heat-inactivate sample at 56°C for 30 minutes in a water bath.
- Allow sample to cool down before retightening cap.
- Sample can be stored frozen until analysis.

Repeated freeze-thawing of sample is not recommended.

#### ADDITIONAL MATERIALS REQUIRED BUT NOT PROVIDED

- Deionized or distilled water
- Disposable gloves
- Rocking platform (designed with a rocking speed range of 12 to 16 cycles per minute, and which moves through a 5° to 10° tilt to wash membranes evenly)
- Pipettors and tips of appropriate volume
- Aspirator with sodium hypochlorite trap
- 56°C water bath (optional)
- Sodium hypochlorite for decontamination

#### PREPARATION OF REAGENTS

- DILUTED WASH BUFFER**
  - DILUTED WASH BUFFER should be **prepared fresh prior to use**.
  - Dilute 1 volume of WASH BUFFER CONCENTRATE (20x) with 19 volumes of reagent grade water. Mix well.
- BLOTTING BUFFER**
  - BLOTTING BUFFER should be **prepared fresh prior to use**.
  - Dilute 1 volume of STOCK BUFFER CONCENTRATE (10x) with 9 volumes of reagent grade water. Mix well.
  - Add 1 g of BLOTTING POWDER to every 20 ml of the diluted STOCK BUFFER prepared in step 2(b) above. Stir to ensure powder dissolves completely.
  - Stir again before dispensing.
- WORKING CONJUGATE SOLUTION**

Note : Prepare solution in polypropylene container / beaker.

  - WORKING CONJUGATE SOLUTION should be **prepared fresh prior to use**.
  - For **RAPID ASSAY PROTOCOL**, prepare WORKING CONJUGATE SOLUTION by diluting CONJUGATE at 1:500 into BLOTTING BUFFER, for example, 10 µl CONJUGATE to 5ml BLOTTING BUFFER.
  - For **OVERNIGHT ASSAY PROTOCOL**, prepare WORKING CONJUGATE SOLUTION by diluting CONJUGATE at 1:1000 into BLOTTING BUFFER, for example, 5 µl CONJUGATE to 5ml BLOTTING BUFFER.
- SUBSTRATE SOLUTION (ready to use)**
  - Dispense directly the required volume from the bottle. Use a clean pipette. Cap tightly after use.

#### AMOUNT OF REAGENTS REQUIRED FOR VARIOUS NUMBER OF STRIPS

Reagents	NUMBER OF STRIPS TO BE USED							
	3	6	9	15	20	27	36	
Diluted Wash Buffer (ml)	60	100	140	240	300	400	600	
Blotting Buffer (ml)	20	40	60	80	100	120	160	
Blotting Powder (g)	1	2	3	4	5	6	8	
Working Conjugate Solution (ml)	7	13	19	31	41	55	73	
Conjugate (µl), Rapid Assay	14	26	38	62	82	110	146	
Conjugate (µl), Overnight Assay	7	13	19	31	41	55	73	
Substrate (ml)	7	13	19	31	41	55	73	

#### ASSAY PROCEDURE - RAPID ASSAY

- Note:** a) Users can use either the rapid or overnight assay to run the tests. HIV bands are more developed and more bands may appear with the overnight assay, but the overall performance of the two assays is the same.
- b) Aspirate all used chemicals and reagents into a trap containing Sodium hypochlorite.
- c) All incubations are to be carried out on a rocking platform.

**Caution:** Some samples cause dark patches on the spot of the strip where they are added. To avoid this problem, one should ensure the following:-

- Sample should be added only after BLOTTING BUFFER is added.
- Tilt the tray slightly by elevating either the top or bottom end of the tray. The Blotting Buffer will flow to the lower end of the tray. Add the sample where the Blotting Buffer is collected. When all the samples are added, return the tray back to its original flat position. Always ensure that the strips are kept wet during the process.
- Alternatively, if tilting the tray is not desired, the samples may be added to the top or bottom end of the well. This way if dark patches showed, the reading of the strip results will not be affected.

#### PROCEDURE:

- Add 2 ml of DILUTED WASH BUFFER to each well. **2 ml**
- Using forceps, carefully remove required number of STRIPS from the tube and place numbered side up into each well. Include strips for Strong Reactive, Weak Reactive and Non-Reactive controls. **2 minutes**
- Incubate the strips for 1 to 2 minutes at room temperature (25 ± 3°C) on a rocking platform (speed of 12 to 16 cycles per minute). Remove buffer by aspiration. (Note: Do not allow the strips to dry. Failure may result in water marks on developed strips for some specimens.) **2 ml**
- Add 2 ml of BLOTTING BUFFER to each well. **20 µl**
- Add 20 µl each of patients' sera or controls to appropriate wells. Care should be taken to ensure specimens are not added directly on the strips. **20 µl**
- Cover the tray with the cover provided and incubate for 1 hour at room temperature (25 ± 3°C) on the rocking platform. **60 minutes**

- Carefully uncover the tray to avoid splashing or mixing of samples. Tilt the tray to aspirate the mixture from the wells. Change aspirator tips between samples to avoid cross-contamination. **3 x 2 ml**

- Wash each strip 3 times with 2 ml of DILUTED WASH BUFFER allowing 5 minutes soak on the rocking platform between each wash. **2 ml**

- Add 2 ml of WORKING CONJUGATE SOLUTION to each well. **60 minutes**

- Cover tray and incubate for 1 hour at room temperature (25 ± 3°C) on the rocking platform. **3 x 2 ml**

- Aspirate CONJUGATE from the wells. Wash as in step 8. **2 ml**

- Add 2 ml of SUBSTRATE SOLUTION to each well. **15 minutes**

- Cover tray and incubate for 15 minutes on the rocking platform. (Note: The reaction can be stopped before 15 minutes if all the bands are visible.) **3 x 2 ml**

- Aspirate the SUBSTRATE and rinse the strips at least three times with reagent grade water to stop the reaction (A dark background can result if washing is insufficient at this step).

- Using forceps, gently remove strips onto paper towels. Cover with paper towels and dry. Alternatively, allow strips to dry in the wells of the tray.

- Mount strips on worksheet (non-absorbent white paper). Do not apply adhesive tape over the developed bands. Observe the bands (See Interpretation of Results) and grade the results. For storage, keep the strips in the dark. **2 ml**

#### ALTERNATIVE PROCEDURE - OVERNIGHT ASSAY

#### PROCEDURE:

- Add 2 ml of DILUTED WASH BUFFER to each well. **2 ml**
- Using forceps, carefully remove required number of STRIPS from the tube and place numbered side up into each well. Include strips for Strong Reactive, Weak Reactive and Non-Reactive controls. **2 minutes**
- Incubate the strips for 1 to 2 minutes at room temperature (25 ± 3°C) on a rocking platform (speed of 12 to 16 cycles per minute). Remove buffer by aspiration. (Note: Do not allow the strips to dry. Failure may result in water marks on developed strips for some specimens.) **2 ml**
- Add 2 ml of BLOTTING BUFFER to each well. **20 µl**
- Add 20 µl each of patients' sera or controls to appropriate wells. **20 µl**
- Cover the tray with the cover provided and incubate overnight (16 - 20 hours) at room temperature (25 ± 3°C) on the rocking platform. **overnight**
- Carefully uncover the tray to avoid splashing or mixing of samples. Tilt the tray to aspirate the mixture from the wells. Change aspirator tips between samples to avoid cross-contamination. **3 x 2 ml**
- Wash each strip 3 times with 2ml of DILUTED WASH BUFFER allowing 5 minutes soak on the rocking platform between each wash.

- Add 2 ml of WORKING CONJUGATE SOLUTION to each well. **2 ml**

- Cover tray and incubate for 30 minutes at room temperature (25 ± 3°C) on the rocking platform. **30 minutes**

- Aspirate CONJUGATE from the wells. Wash as in step 8. **3 x 2 ml**

- Add 2 ml of SUBSTRATE SOLUTION to each well. **2 ml**

- Cover tray and incubate for 15 minutes on the rocking platform. **15 minutes**

- Cover tray and incubate for 1 hour at room temperature (25 ± 3°C) on the rocking platform. **3 x 2 ml**

- Aspirate the SUBSTRATE and rinse the strips at least three times with reagent grade water to stop the reaction (A dark background can result if washing is insufficient at this step).

- Using forceps, gently remove strips onto paper towels. Cover with paper towels and dry. Alternatively, allow strips to dry in the wells of the tray.

- Mount strips on worksheet (non-absorbent white paper). Do not apply adhesive tape over the developed bands. Observe the bands (See Interpretation of Results) and grade the results. For storage, keep the strips in the dark.

SUMMARY OF ASSAY PROTOCOLS			
Reagents	Qty	Room Temp Rapid Assay	Room Temp Overnight Assay
Nitrocellulose strip	1	-	-
Wash Buffer	2 ml	1-2 mins	1-2 mins
Blotting Buffer	2 ml	-	-
Specimen	20 µl	60 mins	Overnight (16 - 20 hours)
Wash Buffer	3 x 2 ml	3 x 5 mins	3 x 5 mins
Conjugate	2 ml	60 mins	30 mins
Wash Buffer	3 x 2 ml	3 x 5 mins	3 x 5 mins
Substrate (Ready to use)	2 ml	15 mins (or less)	15 mins (or less)
Distilled Water	3 x 2 ml	-	-

#### QUALITY CONTROL

We recommend that the Non-Reactive, Strong Reactive and Weak Reactive controls be run with every assay regardless of the number of samples tested. In order for the results obtained from any assay to be considered valid, the following conditions must be met:

- NON-REACTIVE CONTROL**  
No HIV-1 and HIV-2 specific bands should be observed on the Non-Reactive control strips. The band for the serum control should be visible (Fig 1c).
- STRONG REACTIVE CONTROL**  
All relevant molecular weight bands must be evident. Figure 1a provides a guide to the relative positioning of bands visualized with the **MP Diagnostics HIV BLOT 2.2** and permits identification of bands observed for the STRONG REACTIVE CONTROL. The bands are p17, p24, p31, gp41, p51, p55, p66, gp120/gp160. Other bands associated with core antigens (p39, p42) may also be visible. Be careful not to misinterpret these as gp41. The envelope antigens, gp41, gp120/gp160 appear as diffuse bands as they are typical of glycoproteins; p55 viral band may appear faintly on the actual

<b>Precautionary Statements:</b>	P260 Do not breathe dust/fume/gas/mist/vapours/spray. P501 Dispose of contents/ container in accordance with local/regional/national/international regulations.
<b>Supplemental Statements:</b>	EUH210 Safety Data Sheet is available on request
<b>Contains:</b>	0.1% Thimerosal

- Avoid Microbial contamination of reagents when opening and removing aliquots from the original vials or bottles.
- Do not pipette by mouth.
- Handle test specimens, nitrocellulose strips, Reactive, Weak Reactive and Non-Reactive Controls as potentially infectious agents.
- Wear laboratory coats and disposable gloves while performing the assay. Discard gloves in bio-hazard waste-bags. Wash hands thoroughly afterwards.
- It is highly recommended that this assay be performed in a biohazard cabinet.
- Keep materials away from food and drink.
- In case of accident or contact with eyes, rinse immediately with plenty of water and seek medical advice.
- Consult a physician immediately in the event that contaminated materials are ingested or come in contact with open lacerations, or other breaks in the skin.
- Wipe spills of potentially infectious materials immediately with absorbent paper and swab the contaminated area with 1% sodium hypochlorite solution before work is resumed. Sodium hypochlorite should not be used on acid containing spills unless the area is wiped dry with absorbent paper first. Material used (including disposable gloves) should be disposed off as potentially biohazardous material. Do not autoclave material containing sodium hypochlorite.
- Autoclave all used and contaminated materials at 121°C at 15 p.s.i. for 30 minutes before disposal. Alternatively, decontaminate materials in 5% sodium hypochlorite solution for 30-60 minutes before disposal in biohazard waste-bags.
- Decontaminate all used chemicals and reagents by adding sufficient volume of sodium hypochlorite to make a final concentration of at least 1%. Leave for 30 minutes to ensure effective decontamination.
- We do not recommend re-use of incubation trays.

#### ANALYTICAL PRECAUTIONS

- Optimal assay performance requires **STRICT ADHERENCE** to the assay procedure described in this Instructions For Use. Deviations from the procedure may lead to aberrant results.
- DO NOT MODIFY OR SUBSTITUTE REAGENTS FROM ONE KIT LOT TO ANOTHER.** Controls, conjugate and Western Blot strips are matched for optimal performance. Use only the reagents supplied with the kit.
- Do not use kit components beyond the expiry date printed on the kit box.
- Avoid microbial contamination of the reagents when opening and removing aliquots from the original vials or bottles, as this will prematurely reduce the shelf life of the kits and give erroneous results. Use aseptic techniques including pipettes or disposable pipette tips when drawing aliquots from vials.
- The kit controls should be assayed concurrently with patients' samples for each test run.
- Use a new pipette tip for each specimen aliquot to prevent cross contamination.

#### INTERPRETATION OF RESULTS

**NOTE:** Developed strips must be completely dry to avoid misinterpretation.

The presence or absence of antibodies to HIV-1 sample is determined by comparing each nitrocellulose strip to the assay control strips tested with the NON-REACTIVE, STRONG REACTIVE and WEAK REACTIVE controls.

Figure 1a is suggested as an aid to identify the various bands developed on the STRONG REACTIVE Control strip. The Strong Reactive Control as provided in the kit may contain relatively low titer of anti-p55 and anti-p39; as a result, p55 and p39 band for the Strong Reactive Control may appear faintly on the assayed strips. This has no impact on the performance of HIV Blot 2.2 strips in detecting anti-p55 and anti-p39 present in the specimens, as each lot of strip contains sufficient amount of p55 and p39 antigens.

#### PLEASE NOTE:

The numbered end of the strips should be placed at the bottom as shown in the Figure, i.e. the gp120/gp160 bands are the furthest away from the numbered end.

MOLECULAR WEIGHT	GENE	ANTIGEN	DESCRIPTION
gp 160	ENV	Polymeric form of gp41	Broad diffuse glycoprotein
gp 120	ENV	Outer membrane	Diffuse glycoprotein
p66	POL	Reverse Transcriptase	Discreet band
p55	GAG	Precursor protein	Discreet band
p51	POL	Reverse Transcriptase	Discreet band just below p55
p39	GAG	Fragment of p55	Discreet band
gp41	ENV	Transmembrane	Diffuse glycoprotein
p31	POL	Endonuclease	Doublet
p24	GAG	Core protein	Broad band
p17	GAG	Core protein	

PATTERN	INTERPRETATION
No viral specific bands present	NEGATIVE
Detection of p17 antibodies <b>ONLY</b> , no other bands	NEGATIVE
Detection of 2 ENV (gp160/gp41 and gp120) and GAG (p17, p24, p55) or POL (p31, p51, p66)	HIV-1 POSITIVE
Detection of 2 ENV (gp160/gp41 and gp120) and GAG (p17, p24, p55) or POL (p31, p51, p66) and HIV-2 specific band is visible	HIV-1 POSITIVE with HIV-2 INDICATED
Any viral specific bands present but pattern <b>does not</b> meet criteria for POSITIVE	INDETERMINATE <sup>2</sup>
Any viral specific bands present but pattern <b>does not</b> meet criteria for POSITIVE but HIV-2 specific band is visible.	INDETERMINATE <sup>2</sup> with HIV-2 INDICATED

#### INTERPRETATION OF RESULTS FOR INDETERMINATE

INDETERMINATE results should not be used as the basis for diagnosis of HIV-1 infection. Based on the fact that most persons with an initial INDETERMINATE result who are infected with HIV-1 will develop detectable HIV antibodies within 1 month, US CDC (2001) recommended such persons be re-tested for HIV-1 infection  $\geq$  1 month later. Persons with continued INDETERMINATE results after 1 month are unlikely to be HIV-infected unless recent HIV exposure is suspected.

Based on a recent study of Fiebig *et al* (2003), although the window period for Western Blot in the case of a primary HIV-1 infection could be as long as 22 days, the progression from an INDETERMINATE blot to a full POSITIVE profile took no longer than 8 days. In addition, this laboratory stage of having Western Blot INDETERMINATE was always accompanied with detectable RNA of HIV-1 with cases of true infection. Conversely, no seroconversion was evident in follow-up studies of individuals having screened positive and Western Blot INDETERMINATE results, once confirmed as negative by PCR methods (Sethoe *et al*, 1995). Therefore, it is reasonable to consider persons having Western Blot INDETERMINATE results but additionally tested negative by a RNA test as unlikely to be HIV-infected, especially when the tested individuals are known as not having any risk factor associated with exposure.

In particular, persons having Western Blot INDETERMINATE results derived from a test algorithm using fourth generation ELISAs as the primary screen test should additionally be tested for viral RNA using a molecular-base test such as RT-PCR with primer sets covering HIV-1/2/O. If necessary, a follow-up should be considered with any supplemental test 1 month later. The unique design of fourth generation ELISAs is for a simultaneous detection of both antigen and antibody. Consequently specimens identified as positive by a fourth generation ELISA should contain either antibody or antigen or both. Although more than 95% of those cases of true positive identified by a fourth generation ELISA were anti-HIV related and verifiable (confirmed) by Western Blot (Ly *et al.*, 2000), a supplemental test using RT-PCR appeared unavoidable for the small portion of reactivity relating to p24 antigen. Again, persons without any risk of exposure are unlikely HIV-infected, if identified as positive by a fourth generation ELISA accompanied by a Western Blot INDETERMINATE but the findings could not be further supported by a POSITIVE result using a RNA test with primer sets covering HIV-1/2/O.

However, nucleic acid tests (NAT) for HIV DNA or RNA were not approved for diagnostic purpose by the relevant authorities (US CDC, 2001; Constantine & Zink, 2005) until very recently. To date, only one RNA qualitative assay has been approved by the US FDA for diagnosis of primary and acute infection of HIV-1. Therefore, test algorithms recommended by the US CDC (2001) and WHO (2004) are yet to be updated, and NAT are yet to be included as methods for resolving INDETERMINATE Western Blot results. Nevertheless, US CDC (2001) acknowledged that when in consultation with clinical and infection status among persons with an initial INDETERMINATE Western Blot.

#### LIMITATION OF THE METHOD

Detection of antibodies to HIV-1 does not constitute a diagnosis of Acquired Immune Deficiency Syndrome (AIDS). A NEGATIVE BLOT is not a guarantee that the causative agent for AIDS is not present. Although a blot POSITIVE for antibodies to HIV-1 indicates infection with the virus, a diagnosis of AIDS can only be made clinically if a person meets the case definition of AIDS established by the Center for Disease Control (USA), the World Health Organization or other relevant authorities.

It is known that persons who have recently seroconverted may display incomplete pattern but increase reactivity (both number and intensity of bands) occurs when followed for a period of two to six months. Most blots with POSITIVE results will have other viral specific bands present.

INDETERMINATE results should not be used as the basis for diagnosis of HIV-1 infection. It is recommended that all INDETERMINATE blots be repeated using the original specimen and sequential samples. Blood donors with an INDETERMINATE blot should be re-tested using a fresh specimen after one month (US CDC, 2001). In addition, antibodies to p24 and p31 are known to decrease during the course of AIDS leading to a shift in blot interpretation from POSITIVE to INDETERMINATE. Interpretation of results should then be based on subsequent blot testing and clinical evaluations in such situations.

Due to its highly specific nature, NON-REACTIVITY of samples with HIV-2 specific envelope peptide on an Indeterminate viral blot, does not exclude the possibility of infection with other strains of HIV-2.

Samples that are indicated as HIV-2 infections should be further tested with a HIV-2 Western Blot Kit.

#### SPECIFIC PERFORMANCE CHARACTERISTICS

The performance of MP Diagnostics HIV BLOT 2.2 for the detection of antibodies to HIV-1 or HIV-2 was evaluated in clinical studies.

Table 1 : Sensitivity study of HIV-1 viral antigen reactivity with HIV-1 seropositive samples. (Number of samples = 201)

SEROLOGICAL PROFILE	HIV BLOT 2.2	DUPONT/ORTHO HIV-1 WB
GAG, POL and ENV	97.5%	95.4%
p24, p31, gp41 and/or gp120/gp160	94.9%	90.9%
ENV and GAG or POL	100.0%	100.0%

Table 2: Specificity study of HIV-1 viral antigen reactivity with normal donor samples and sera with other viral infections.

SAMPLE TYPE	NUMBER TESTED	POSITIVE	HIV-1 REACTIVITY	
			INDETERMINATE <sup>1</sup>	NEGATIVE
Normal Donors	208	0	11	197
HTLV-1	5	0	0	5
CMV	5	0	1	4
EBV (IgM)	5	0	1	4
V.zoster (IgG)	5	0	1	4
Measles	6	0	2	4
Rubella	5	0	1	4
Mumps	4	0	1	3
Adenovirus	5	0	2	3
HSV	5	0	0	5
Dengue	5	0	1	4
<b>Total</b>	<b>258</b>	<b>0</b>	<b>21</b>	<b>237</b>

<sup>1</sup>All showed as a p24 or p17 band only.

Table 3 : Sensitivity study of HIV-2 peptide band with HIV-2 seropositive samples. (Number of samples = 178)

HIV-2 Western Blot Serological profile*	HIV-2 peptide Reactivity	
	Positive	Negative
GAG, POL and 2 ENV	160	0
GAG, POL and 1 ENV	18	0

\*Sera define as positive by results of Pasteur New LAV Blot 2. Data provided by Dr. Oliviero E. Varnier and Dr. Flavia Lillo. Laboratory of Human Retroviruses. University of Genoa.

Table 4 : Specificity study of HIV-2 peptide band with HIV-1 seropositive sera, normal donor samples and sera with other viral infections.

SAMPLE TYPE	NUMBER TESTED	HIV-2 PEPTIDE REACTIVITY	
		POSITIVE	NEGATIVE
HIV-1 seropositive	197	16*	181
Normal Donors	208	0	208
HTLV-1 seropositive	5	0	5
CMV	5	0	5
EBV (IgM)	5	0	5
V.zoster (IgG)	5	0	5
Measles	6	0	6
Rubella	5	0	5
Mumps	4	0	4
Adenovirus	5	0	5
HSV	5	0	5
Dengue	5	0	5
<b>Total</b>	<b>455</b>	<b>16</b>	<b>439</b>

\*When tested on the MP Diagnostics HIV-2 Western Blot, 6 of these samples had reactivity with ENV and GAG or POL, and 9 of these samples had reactivity to only GAG and/or POL while 1 sample was negative.

A total of 15 commercial HIV-1 seroconversion panels were tested with MP Diagnostics HIV Blot 2.2 and results showed that the MP Diagnostics HIV Blot 2.2 was able to detect antibody to HIV earlier or in the same sample in all the panels.

#### LIMITED EXPRESSED WARRANTY DISCLAIMER

The manufacturer makes no expressed warranty other than that the test kit will function as an *in vitro* diagnostic assay within the specifications and limitations described in the Product Instructions For Use when used in accordance with the instructions contained therein. The manufacturer disclaims any warranty, expressed or implied, including such expressed or implied warranty with respect to merchantability, fitness for use or implied utility for any purpose. The manufacturer is limited to either replacement of the product or refund of the purchase price of the product. The manufacturer shall not be liable to the purchaser or third parties for any damage, injury or economic loss howsoever caused by the product in the use or in the application thereof.

#### TECHNICAL PROBLEMS / COMPLAINTS

Should there be a technical problem / complaint, please do the following :

- Note the kit lot number, the expiry date and the strip lot number.
- Retain the kits and the results that were obtained.
- Contact the nearest MP Biomedicals office or your local distributor.

#### BIBLIOGRAPHY

- V.C.W.Tsang, K. Hancock, M. Wilson. D.F. Palmer, S. Whaley, J.S. Mc Dougal, and S. Kennedy. March 1985. Developmental Procedure : Enzyme-linked Immuno-electro-transfer Blot technique for HTLV-III/LAV antibodies; CDC, Atlanta.
- H. Towbin, T. Staehlin, and J. Gordon. 1979. Electrophoretic transfer of proteins from polyacrylamide gels to nitrocellulose sheets: procedure and some applications. Proc. Natl. Acad. Sci. USA 76: 4350-4354.
- J. Schupbach, M. Popovic, R. V. Gilden. M.A. Gonda, M. G. Sarnagadharan and R. C. Gallo. 1984. Serological Analysis of subgroup of Human T-Lymphotropic retroviruses (HTLV-III) associated with AIDS. Science 224, 503-505.
- M. G. Sarnagadharan, M. Popovic, L. Bruch, J. Schupbach and R. C. Gallo. 1984. Antibodies reactive with human T-Lymphotropic retroviruses (HTLV-III) in the serum of patients with AIDS. Science 224, 506-608.
- CDC. 1985. " Provisional public health service inter-agency recommendations for screening donated blood and plasma for antibody to the virus causing Acquired Immune Deficiency Syndrome" - United States Morbidity and Mortality Weekly Report 34 (1) :1-5.
- Proposed World Health Organization 1990 criteria for interpreting results from Western blot assays for HIV-1, HIV-2, and HTLV-I/HTLV-II, Weekly Epidemiological Record 65(37), 281-283.
- F. Clavel, D. Guetard., F. Brun-Vezinet, et al. 1986 Isolation of a new human retrovirus from West African patients with AIDS. Science; 233:343-346.
- F. Clavel., 1987. HIV-2, the West African AIDS virus. AIDS 1:135-140.
- R.S. Tedder, A. Hughes, T. Corrah et al 1988. Envelope cross-reactivity in Western Blot for HIV-1 and HIV-2 may not indicate dual infection. Lancet 11:927-930.
- Bottiger B., A. Karlsson, F. Andreasson et al. 1990. Envelope cross-reactivity between Human Immunodeficiency Virus Type 1 and Type 2 detected by different serological methods: Correlation between cross-neutralization and reactivity against the main neutralizing site. J. Virol. 64(7):3492-3499.
- Centers for Disease Control. 2001. Revised Guidelines for HIV Counseling, Testing, and Referral and Revised Recommendations for HIV Screening of Pregnant Women – United States, Morbid. Mortal. Weekly Rep. 50: RR-19.
- Fiebig, E. W., D. J. Wright, B. D. Rawal, P. E. Garrett, R. T. Schumacher, L. Peddada, C. Heldebrandt, R. Smith, A. Conrad, S. H. Kleinman, and M. P. Busch. 2003. Dynamics of HIV viremia and antibody seroconversion in plasma donors: implications for diagnosis and staging of primary HIV infection. AIDS. 17:1871-1879.
- Ly, T. D., C. Edlinger, and A. Vabret. 2000. Contribution of combined detection assays of p24 antigen and anti-human immunodeficiency virus (HIV) antibodies in diagnosis of primary HIV infection by routine testing. J Clin Microbiol. 38:2459-2461.
- Sethoe, S. Y., A. E. Ling, E. H. Sng, E. H. Monteiro, R. K. Chan. 1995. PCR as a confirmatory test for human immunodeficiency virus type 1 infection in individuals with indeterminate western blot (immunoblot) profiles. J Clin Microbiol. 33:3034-3036.
- Constantine, N. T. and H. Zink. 2005. HIV testing technologies after two decades of evolution. Indian J Med Res. 121:519-538.
- World Health Organization. 2004. Guidelines for HIV Diagnosis and monitoring of antiretroviral therapy. Regional Office for South-East Asia, New Delhi, India.
- Ming Guan, Frequency, causes and new challenges of indeterminate results in Western Blot Confirmatory Testing for Antibodies to Human Immunodeficiency Virus. Clinical and Vaccine Immunology, June 2007, Vol.14, No.6, p649-659.

**MP Biomedicals Asia Pacific Pte Ltd.**  
 2 Pioneer Place  
 Singapore 627885  
 Tel. No. : + 65 6775 0008  
 Fax. No. : + 65 6774 6146  
 Email : enquiry\_ap@mpbio.com


**EC REP** **MP Biomedicals Germany GmbH**  
 Thüringer Straße 15  
 37269 Eschwege  
 Germany  
 Tel. No. : +49 5651 921 204  
 Fax No. : +49 5651 921 181  
 Email : diagnostics@mpbio.com

#### Regional Office:

**MP Biomedicals Germany GmbH**  
 Thüringer Straße 15  
 37269 Eschwege  
 Germany  
 Tel. No. : +49 5651 921 204  
 Fax No. : +49 5651 921 181  
 Email : diagnostics@mpbio.com

\* U.S. Patent 5,721,095

FIGURE 1



- Strong Reactive Control (Reactive for HIV-1 and HIV-2)
- Weak Reactive Control (Reactive for HIV-1 only).
- Non-Reactive Control.
- A typical HIV-2 seropositive serum.

#### TROUBLE SHOOTING CHART

