

# MP Biomedicals

## Multitest Slides

One Call

Build welcoming homes  
for your cells!

One Source

A World of  
Cell Culture



Many applications including immunofluorescence,  
cell culture, slide agglutination tests and histology

These high quality glass slides have a unique PTFE coating which offers several important advantages:

- ▶ Self draining surface
- ▶ Cells only attach to uncoated glass
- ▶ Resistant to common organic solvents and fixatives
- ▶ Resistant to direct heat




[www.mpbio.com](http://www.mpbio.com)

## Description

The slides are available in 8, 10, 12 and 15 well format while the 8 well slides feature the same spacing as standard microplates. The slide size is 2.5 cm x 7.5 cm x 1 mm (100 /Box). Easily sterilized in 1 hour at 160 °C.

Cat. No	Pack Size	No. of well	Well Size
096040805E	100/box	8	6 mm
096041805E		10	8 mm
096041205E		12	8 mm
096041505E		15	4 mm



## Principle

These high quality glass slides present a blue single PTFE covering, ensuring speed of drying, resistance to principal organic solvents and possibility of sterilization by high temperature.

The uncoated areas have a total wettable surface so that liquid spreads evenly over the area of each well.

For cell attachment and growth, the slides require only a simple washing in deionized water or detergent followed by sterilization.



- Homogeneous distribution of cellular equipment
- Several simultaneous immunofluorescences (IF)
- Multiparametric analysis
- High throughput (compatible multichannel pipette)
- Agglutination tests
- Resistant to organic solvents
- Autoclavables

**Order Multitest Slide for Your Lab Today!**

<b>MP Biomedicals EMEA</b> MP Biomedicals SARL Parc d'Innovation, BP 50067 67402 Illkirch Cedex, France Tel: 33 03 8867 5425 Toll Free: 800 7777 9999 Fax: 33 03 8867 1945 Email: custserv.eur@mpbio.com	<b>MP Biomedicals Worldwide</b> MP Biomedicals LLC, HQ 3 Hutton Center Drive Suite 100, Santa Ana, California 92707 Tel: 1 949 833 2500 Toll Free: 1 800 854 0530 Fax: 1 949 859 5010 Email: sales@mpbio.com	<b>MP Biomedicals APAC</b> MP Biomedicals Asia Pacific Pte Ltd, HQ 2 Pioneer Place Singapore 627885 Tel: 65 6775 0008 Fax: 65 6774 6146 Email: enquiry_ap@mpbio.com	<b>For other offices and more information, please visit <a href="http://www.mpbio.com">www.mpbio.com</a></b>
---	---	--	--