



MP Biomedicals, LLC

29525 Fountain Parkway
Solon, Ohio 44139

Telephone: 440/337-1200
Toll Free: 800/854-0530
Fax: 440/337-1180
mailto: biotech@mpbio.com
web: <http://www.mpbio.com>

TECHNICAL INFORMATION

Catalog Number: 191391

Complement C1Q Protein

CAS # 80295-33-6

Source: Human Serum

Physical Description: Supplied in PBS (10 mM sodium phosphate, 135 mM sodium chloride, pH ~7.2, 40% glycerol).

Concentration: ~ 1 mg/ml

Recommended Storage: -70°C

Note: All raw materials derived from human origin used in the manufacture of this product have been tested for and found to be non-reactive for HbsAG, HCV, HIV I and II by currently approved FDA methods. We certify this to be true to the best of our knowledge. However, because no test method can offer complete assurance that HIV, HCV or other infectious agents are absent, this material should be handled as a biosafety hazard.