



MP Biomedicals, LLC

29525 Fountain Parkway  
Solon, Ohio 44139

Telephone: 440/337-1200  
Toll Free: 800/854-0530  
Fax: 440/337-1180  
mailto: [biotech@mpbio.com](mailto:biotech@mpbio.com)  
web: <http://www.mpbio.com>

## TECHNICAL INFORMATION

# ImmunO™

Catalog Number: 643711

### Rabbit Anti-Dog IgM Serum (Lyophilized)

**Form:** Lyophilized antiserum.

**Description:** Prepared in rabbits with highly purified normal dog IgM in Freund's adjuvant. Then the antiserum is solid phase absorbed to ensure specificity to the mu chain.

**Reconstitution:** Reconstitute with 2.0 ml of distilled water.

**Specificity:** This product is monospecific for dog IgM by IEP and gel diffusion analysis.

**Note:** This product may contain a preservative such as sodium azide, thimerosal or proclin. Please see lot specific chemical credential for preservative information.

[If a titer/working dilution is not given above, please click here to see a general dilution chart for working with antibodies. Please note that the general dilution chart should only be used as a guideline. Each lab should determine their own optimal working dilution.](#)

[Will this antibody work with your application? Please click here to see a general chart of antibody applications. Please note that any information given above is primary application data. The general applications charts should only be used as a reference.](#)