



MP Biomedicals, LLC

29525 Fountain Parkway  
Solon, Ohio 44139

Telephone: 440/337-1200  
Toll Free: 800/854-0530  
Fax: 440/337-1180  
mailto: [biotech@mpbio.com](mailto:biotech@mpbio.com)  
web: <http://www.mpbio.com>

## TECHNICAL INFORMATION

# cappel

Catalog Number: 56963

### Rabbit Affinity Purified Antibody To Bovine Albumin

**DESCRIPTION:** Product is rabbit affinity purified antibody to bovine albumin and buffer salts.

The affinity purified antibody is prepared from the specific rabbit IgG fraction by solid phase absorption. The resulting affinity purified antibody is dialyzed into 0.02M sodium phosphate, 0.14M sodium chloride, pH 7.3, adjusted to standard titer, filtered through a 0.22  $\mu$ m filter, and vialled.

**USES:** This product is suitable for use as a precipitating antibody. It gives sharp precipitin lines in immunoelectrophoresis and immunodiffusion techniques.

This product is suitable for conjugation with enzymes, radiolabels, fluorochromes or other markers, and for attachment to solid supports for use as an immunoabsorbant.

This product is suitable for use as a primary reagent in enzyme immunoassays, fluorescein immunoassays and blot, cell or tissue immunostaining. These affinity purified antibodies are recommended to avoid non-specific binding from the inherent antibodies of host animals.

See **Related Product Ordering Information**.

**SPECIFICATIONS:** The total protein is measured by absorbance ( $A_{280nm}$ ). Antibody titer is standardized by equivalence-point precipitation. Each affinity purified antibody is tested for purity at 10 mg/ml using immunoelectrophoresis. The product is rabbit IgG; no trace of albumin is detected. It shows reactivity to bovine albumin; cross-reactivity to other species may exist. Antibody activity to other serum proteins is not present.

#### RELATED PRODUCT ORDERING INFORMATION:

Catalog #DescriptionSize

55269 RABBIT ANTI-BOVINE ALBUMIN2ML  
55275 IGG RABBIT ANTI-BOVINE ALBUMIN 5ML  
55281 FITC RABBIT ANTI-BOVINE ALBUMIN 2ML  
55283 RHO RABBIT ANTI-BOVINE ALBUMIN2ML  
55285 HRP RABBIT ANTI-BOVINE ALBUMIN2ML  
55886 FITC BOVINE ALBUMIN 25MG  
55897 RHODAMINE BOVINE ALBUMIN 25MG  
55918 PURIFIED BOVINE ALBUMIN50MG

**Note:** This product may contain a preservative such as sodium azide, thimerosal or proclin. Please see lot specific chemical credential for preservative information.

[If a titer/working dilution is not given above, please click here to see a general dilution chart for working with antibodies. Please note that the general dilution chart should only be used as a guideline. Each lab should determine their own optimal working dilution.](#)

[Will this antibody work with your application? Please click here to see a general chart of antibody applications. Please note that any information given above is primary application data. The general applications charts should only be used as a reference.](#)